pa_wkit
water activity meter

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8. Odkazy na literaturu

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- Bousquet-Ricard, M., G. Qualyle, T. Pharm, and J. C. Cheftel. 1980. Comparative study of three methods of determining water activity in intermediate moisture foods. *Lebensm Wiss Technol* 13:169-173.
- Cazier, J.B., and V. Gekas. 2001. Water activity and its prediction: a review. *International Journal of Food properties* 4(1):35-43.
- Chirife, J., G. Favetto, C. Ferro-Fontn, and S.L.Resnik. 1983. The water activity of standard saturated salt solutions in the range of intermediate moisture foods. *Lebensm Wiss Technol* 16:36-38.
- Duckworth, R. (1975). *Water Relations of Foods*. Academic Press, New York.
- Gomez-Diaz, R. (1992). Water activity in foods: Determination methods. *Alimentaria*. 29:77-82.
- Gmez, R., and J. Fernandez-Salguero. 1992. Water activity and chemical composition of some food emulsions. *Food Chem* 45:91-93.
- Greenspan, L. (1977). Humidity fixed points of binary saturated aqueous solutions. *Journal of Research of the National Bureau of Standards - A.Physics and Chemistry*. 81A:89-96.
- Karmas, E. 1981. Measurement of moisture content. *Cereal Foods World* 26:332-334.
- Kitic, D., D.C. Pereira-Jardim, G.J. Favetto, S.L. Resnik, and J. Chirife. 1986. Theoretical prediction of the water activity of standard saturated salt solutions at various temperatures. *Journal of Food Science* 51:1037-1042.
- Labuza, T.P., and R. Contreras-Medellin. 1981. Prediction of moisture protection requirements for foods. *Cereal Foods World* 26:335-343.
- Labuza, T.P., K. Acott, S.R. Tatini, R.Y. Lee, J. Flink, and W. McCall. 1976. Water activity determination: A collaborative study of different methods. *Journal of Food Science* 41:910-917.
- Marculli, C., and Th. Peter. 2005. Water activity in polyol/water systems: new UNIFAC parameterization. *Atmospheric Chemistry and Physics* 5:1545-1555.
- Ninni, L., M.S. Camargo, and A.J.A. Meirelles. 2000. Water activity in polyol systems. *Journal of Chemical and Engineering Data* 45:654-660.
- Prior, B.A. (1979). Measurement of water activity in foods: A review. *Journal of Food Protection*. 42(8):668-674.
- Rahman, M.S. and S.S. Sablani. 2001. Measurement of water activity by electronic sensors. P. A2.5.1-A2.5.4 In R.E.Wrolstad (ed.) *Current Protocols In Food Analytical Chemistry*. John Wiley & Sons, Inc., New York.
- Rahman, M.S., S.S. Sablani, N. Guizani, T.P. Labuza, and P.P. Lewicki. 2001. Direct manometric determination of vapor pressure. P. A2.4.1-A2.4.6. In R.E. Wrolstad (ed.) *Current Protocols In Food Analytical Chemistry*. John Wiley & Sons, Inc., New York.
- Reid, D.S., A.J. Fontana, M.S. Rahman, S.S. Sablani, T.P. Labuza, N. Guizani, and P.P. Lewicki. 2001. Vapor pressure measurements of water p. A2.1.1-A2.5.4. In R.E. Wrolstad (ed.) *Current Protocols In Food Analytical Chemistry*. John Wiley & Sons, Inc., New York.

- Reid, D.S. 1976. Water activity concepts in intermediate moisture foods. p. 54-65. In R. Davies, G.G. Birch, and K.J. Parker (ed.) *Intermediate Moisture Foods*. Applied Science Publishers, London.
- Richard, J., and T.P. Labuza. 1990. Rapid determination of the water activity of some reference solutions, culture media and cheese using a dew point method. *Sci. des Aliments* 10:57-64.
- Roa, V., and M.S. Tapia de Daza. 1991. Evaluation of water activity measurements with a dew point electronic humidity meter. *Lebensm Wiss Technol* 24:208-213.
- Rodel, W. 2001. Water activity and its measurement in food. P. 453-483. In E. Kress-Rogers, and C.B. Brimelow (ed.) *Instrumentation and sensors for the food industry*. CRC Press LLC, Boca Raton, FL.
- Roos, K.D. 1975. Estimation of water activity in intermediate moisture foods. *Food Tech* 29:26-30.
- Scott, V.N., and D.T. Bernard. 1983. Influence of temperature on the measurement of water activity of food and salt systems. *Journal of Food Science* 48:552-554.
- Snavely, M.J., J.C. Price, and H.W. Jun. 1990. A comparison of three equilibrium relative humidity measuring devices. *Drug Dev. Ind. Pharm.* 16:1399-1409.
- Stamp, J.A., S. Linscott, C. Lomauro, and T.P. Labuza. 1984. Measurement of water activity of salt solutions and foods by several electronic methods as compared to direct vapor pressure measurement. *Journal of Food Science* 49:1139-1142.
- Stoloff, L. 1978. Calibration of water activity measuring instruments and devices: Collaborative study. *Journal of the Association of Official Analytical Chemists* 61:1166-1178.
- Troller, J.A. 1983. Methods to measure water activity. *Journal of Food Protection* 46:129-134.
- Troller, J.A. and J.H.B. Christian. (1978). *Water Activity and Food*. Academic Press, New York.
- Troller, J.A. and V.N. Scott. (1992). Measurement of water activity (A_w) and acidity. In: *Compendium of Methods for the Microbiological Examination of Foods*. Vanderzant, C. and D.F. Splittstoesser (ed.) American Public Health Association, Washington, D.C. pp. 135-151.
- van den Berg, C. (1985). Water activity. In: *Concentration and Drying of Foods*. MacCarthy, D. (ed.) Elsevier, London. pp. 11-35.
- Van den Berg, C. 1986. Water activity. p. 11-36. In D. MacCarthy (ed.) *Concentration and drying of foods*. Elsevier Applied Science Publishers, London.
- Van den Berg, C. 1991. Food-water relations: Progress and integration, comments and thoughts. In H. Levine, and L. Slade (ed.) *Water Relationships in Foods*. Plenum Press, New York.
- Van den Berg, C., and S. Bruin. 1981. Water activity and its estimation in food systems: Theoretical aspects. p. 1-61. In L.B. Rockland, and G.F. Stewart (ed.) *Water Activity: Influences on Food Quality*. Academic Press, New York.
- Vega-Mercado, H., and G.V. Barbosa-Canovas. 1994. Prediction of water activity in food systems: A review on theoretical models. *Revista Espanola De Ciencia Y Tecnologia De Alimentos* 34:368-388.
- Vega-Mercado, H., B. Romanach, and G.V. Barbosa-Canovas. 1994. Prediction of water activity in food systems: A computer program for predicting water

- activity in multicomponent foods. *Revista Espanola De Ciencia Y Tecnologia De Alimentos* 34:427-440.
- Vos, P.T., and T.P. Labuza. 1974. Technique for measurements of water activity in the high aw range. *J. Agric. Food Chem.* 22:326-327.
- Voysey, P. 1993. An evaluation of the AquaLab CX-2 system for measuring water activity. *F. M. B. R. A. Digest No.* 124 24-25.

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- Bei, Z.H., and R.-M.J. Nout. 2000. Effects of temperature, water activity and gas atmosphere on mycelial growth of tempe fungi *Rhizopus microsporus* var. *microsporus* and *R. microsporus* var. *oligosporus*. *World Journal of Microbiology and Biotechnology* 16:853-858.
- Beuchat, L.R. 1981. Microbial stability as affected by water activity. *Cereal Foods World* 26:345-349.
- Brandt, L. (1996). Bound for success. Controlling water activity gives technologists the edge in developing safe, shelf-stable foods. *Food Formulating*. September:41-48.
- Chirife, J., and M.P. Buera. 1994. Water activity, glass transition and microbial stability in concentrated/semimoist food systems. *Journal of Food Science* 59:921-927.
- Chirife, J., and M.P. Buera. 1995. A critical review of some nonequilibrium situations and glass transitions on water activity values of foods in the microbiological growth range. *Journal of Food Engineering* 25:531-552.
- Chirife, J., and M.P. Buera. 1996. Water activity, water glass dynamics, and the control of microbiological growth in foods. *Critical Rev. in Food Sci. Nutr.* 36:465-513.
- Farber, J.M., F. Coates, and E. Daley. 1992. Minimum water activity requirements for the growth of *Listeria monocytogenes*. *Lett Appl Microbiol* 15:103-105.
- Franks, F. 1991. Water activity: a credible measure of food safety and quality? *Trends Food Sci Technol* March:68-72.
- Franks, F. (1982). Water activity as a measure of biological viability and quality control. *Cereal Foods World*. 27(9):403-407.
- Garcia de Fernando, G.D., O. Diaz, M. Fernandez, and J.A. Ordonez. 1992. Changes in water activity of selected solid culture media throughout incubation. *Food Microbiology* 9:77-82.
- Gibson, A.M., J. Baranyi, J.I. Pitt, M.J. Eyles, and T.A. Roberts. 1994. Predicting fungal growth: The effect of water activity on *Aspergillus avus* and related species. *International Journal of Food Microbiology* 23:419-431.
- Goeleni, N., J.E. Smith, J. Lacey, and G. Gettinby. 1997. Effects of temperature, water activity, and incubation time on production of aatoxins and cyclopiazonic acid by an isolate of *Aspergillus avus* in surface agar culture. *Appl Environ Microbiol* 63:1048-1053.
- Hardman, T.M. (1988). *Water and Food Quality*. Elsevier Press, London.
- Hocking, A.D., B.F. Miscamble, and J.I. Pitt. 1994. Water relations of *Alternaria alternata*, *Cladosporium cladosporioides*, *Cladosporium sphaerospermum*, *Curvulario lunata* and *Curvulario pallescens*. *Mycological Research* 98:91-94.
- Houtsma, P.C., A. Heuvelink, J. Dufrenne, and S. Notermans. 1994. Effect of sodium lactate on toxin production, spore germination and heat resistance

- of proteolytic *Clostridium botulinum* strains. *Journal of Food Protection* 57:327-330.
- Kress-Rogers, E. (1993). Food quality measurement. *Food Industry News*. September:23-26.
- Kuntz, L.A. 1992. Keeping microorganisms in control. *Food Product Design* August:44-51.
- Levine, H., and L. Slade. 1991. *Water Relationships in Foods*. Plenum Press, New York.
- Li, K.Y., and J.A. Torres. 1993. Water activity relationships for selected mesophiles and psychrotrophs at refrigeration temperature, *Journal of Food Protection* 56:612-615.
- Lopez-Malo, A., S. Guerrero, and S.M. Alzamora. 2000. Probabilistic modeling of *Saccharomyces cerevisiae* inhibition under the effects of water activity, pH, and potassium sorbate concentration. *Journal of Food Protection* 63:91-95.
- Mannheim, C.H., J.X. Liu, and S.G. Gilbert. 1994. Control of water in foods during storage. *Journal of Food Engineering* 22:509-532.
- Marauska, M., A. Vigants, A. Klincare, D. Upite, E. Kaminska, and M. Bekers. 1996. Influence of water activity and medium osmolality on the growth and acid production of *Lactobacillus casei* var. *alactosus*. *Proceedings of the Latvian Academy of Sciences Section B. Natural Exact and Applied Sciences* 50:144-146.
- Masana, M.O., and J. Baranyi. 2000. Growth/no growth interface of *Brochothrix thermosphacta* as a function of pH and water activity. *Food Microbiology* 17:485-858.
- Mattick, K. L., F. Jorgensen, J.D. Legan, M.B. Cole, J. Porter, H.M. Lappin-Scott, and T.J. Humphrey. 2000. Survival and filamentation of *Salmonella enterica* serovar *Enteritidis* PT4 and *Salmonella enterica* serovar *Typhimurium* DT104 at low water activity. *Appl Environ Microbiol* 66:1274-1279.
- Mattick, K.L., F. Jorgensen, J.D. Legan, H.M. Lappin-Scott, and T.J. Humphrey. 2000. Habituation of *Salmonella* spp. at reduced water activity and its effect on heat tolerance. *Appl Environ Microbiol* 66:4921-4925.
- Mattick, K.L., F. Jorgensen, J.D. Legan, H.M. Lappin-Scott, and T.J. Humphrey. 2001. Improving recovery of *Salmonella enterica* Serovar *Typhimurium* DT104 cells injured by heating at different water activity values. *Journal of Food Protection* 64:1472-1476.
- McMeekin, T.A. and T. Ross. (1996). Shelf life prediction: Status and future possibilities. *International Journal of Food Microbiology*. 33:65-83.
- Miller, A.J. 1992. Combined water activity and solute effects on growth and survival of *Listeria monocytogenes*. *Journal of Food Protection* 55:414-418.
- Nakajo, M., and Y. Moriyama. 1993. Effect of pH and water activity on heat resistance of spores of *Bacillus coagulans*. *Journal of the Japanese Society for Food Science and Technology* 40:268-271.
- Nelson, K.A., and T.P. Labuza. 1994. Water activity and food polymerscience: Implications of state on arrhenius and WLF models in predicting shelf life. *Journal of Food Engineering* 22:271-289.
- Nesci, A., M. Rodrigues, and M. Etcheverry. 2003. Control of *Aspergillus* growth and aatoxin production using antioxidants at different conditions of water activity and pH. *Journal of Applied Microbiology* 95:279-287.
- Nolan, D.A., D.C. Chamblin, and J.A. Troller. 1992. Minimal water activity levels for growth and survival of *Listeria monocytogenes* and *Listeria innocua*. *International Journal of Food Microbiology* 16:323-335.

- Noorlidah, A., A. Nawawi, and I. Othman. 2000. Fungal spoilage of starch-based foods in relation to its water activity (aw). *Journal of Stored Products Research* 36:47-54.
- Park, C.M., and L.R. Beuchat. 2000. Survival of *Escherichia coli* O157:H7 in potato starch as affected by water activity, pH and temperature. *Lett Appl Microbiol* 31(5):364-367.
- Petersson, S., and J. Schnuerer. 1995. Biocontrol of mold growth in high-moisture wheat stored under airtight conditions by *Pichia anomala*, *Pichia guilliermondii*, and *Saccharomyces cerevisiae*. *Appl Environ Microbiol* 61:1027-1032.
- Pitt, J.I., and B.F. Miscalable. 1995. Water relations of *Aspergillus avus* and closely related species. *Journal of Food Protection* 58:86-90.
- Plaza, P., J. Usall, N. Teixido, and I. Vinas. 2003 Effect of water activity and temperature on germination and growth of *Penicillium digitatum*, *P. italicum* and *Geotrichum candidum*. *Journal of Applied Microbiology* 94:549-554.
- Quintavalla, S., and G. Parolari. 1993. Effects of temperature, water activity and pH on the growth of *Bacillus* cells and spore: A response surface methodology study. *International Journal of Food Microbiology* 19:207-216.
- Rockland, L.B. and G.F. Stewart. (1981). *Water Activity: Influences on Food Quality*. Academic Press, New York.
- Rockland, L.B., and S.K. Nishi. 1980. Influence of water activity on food product quality and stability. *Food Tech* 34:42-59.
- Saad, R.R. 1992. Effect of water activity on growth and lipids of xerophilic fungi, *Aspergillus repens* and *Aspergillus amstelodami*. *Zentralblatt Fuer Mikrobiologie* 147:61-64.
- Salter, M.A., D.A. Ratkowsky, T. Ross, and T.A. McMeekin. 2000. Modelling the combined temperature and salt (NaCl) limits for growth of a pathogenic *Escherichia coli* strain using nonlinear logistic regression. *International Journal of Food Microbiology* 61:159-167.
- Santos, J., T.M. Lopez-Diaz, M.C. Garcia-Lopez, M.C. Garcia-Fernandez, and A. Otero. 1994. Minimum water activity for the growth of *Aeromonas hydrophila* as affected by strain, temperature and humectant. *Lett Appl Microbiol* 19:76-78.
- Sautour, M., A. Rouget, P. Dantigny, C. Divies, and M. Bennisoussan. 2001. Prediction of conidial germination of *Penicillium chrysogenum* as influenced by temperature, water activity and pH. *Lett Appl Microbiol* 32:131-134.
- Seow, C.C., T.T. Teng, and C.H. Quah. (1988). *Food Preservation by Moisture Control*. Elsevier, New York.
- Shebuski, J.R., O. Vilhelmsson, and K.J. Miller. 2000. Effects of growth at low water activity on the thermal tolerance of *Staphylococcus aureus*. *Journal of Food Protection* 63:1277-1281.
- Taoukis, P., W. Breene, and T.P. Labuza. (1988). Intermediate moisture foods. *Advances in Cereal Science and Technology*. 9:91-128.
- Tapia de Daza, M.S., Y. Villegas, and A. Martinez. 1991. Minimal water activity for growth of *Listeria monocytogenes* as affected by solute and temperature. *International Journal of Food Microbiology* 14:333-337.
- Tokuoka, K., and T. Ishitani. 1991. Minimum water activities for the growth of yeasts isolated from high-sugar foods. *Journal of General and Applied Microbiology* 37:111-119.

- Torres, R., J. Usall, N. Teixido, M. Abadias, and I. Vinas. 2003. Liquid formulation of the biocontrol agent *Candida sake* by modifying water activity or adding protectants. *Journal of Applied Microbiology* 94:330-339.
- Ucar, F., and I. Guneri. 1996. The effect of water activity, pH and temperature on the growth of osmophilic yeasts. *Turkish Journal of Biology* 20:37-46.
- Wijtzes, T., P.J. McClure, M.H. Zwietering, and T.A. Roberts. 1993. Modelling bacterial growth of *Listeria monocytogenes* as a function of water activity, pH and temperature. *International Journal of Food Microbiology* 18:139-149.
- Zwietering, M.H., T. Wijtzes, J.C. de Wit, and K. Van't Riet. 1992. A decision support system for prediction of the microbial spoilage in foods. *Journal of Food Protection* 55:973-979.

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8.3.1 Maso a mořské produkty

- Allen, K., D. Cornforth, D. Whittier, M. Vasavada, and B. Nummer. 2007. Evaluation of high humidity and wet marinade methods for pasteurization of jerky. *Journal of Food Science*. 72:C351-C355.
- Chen, H.C. 1995. Seafood microorganisms and seafood safety. *Journal of Food and Drug Analysis* 3:133-144.
- Clavero, M.R.S., and L.R. Beuchat. 1996. Survival of *Escherichia coli* O157:H7 in broth and processed salami as influenced by pH, water activity, and temperature and suitability of media for its recovery. *Appl Environ Microbiol* 62:2735-2740.
- Duffy, L.L., P.B. Vanderlinde, and F.H. Grau. 1994. Growth of *Listeria monocytogenes* on vacuum-packed cooked meats: Effects of pH, aw, nitrite and ascorbate. *International Journal of Food Microbiology* 23:377-390.
- Elgasim, E.A., and M.S. Al Wesali. 2000. Water activity and Hunter colour values of beef patties extended with samh (*Mesembryanthemum forsskaei* Hochst) our. *Food Chem* 69(2):181-185.
- Gmez, R., and J. Fernandez-Salguero. 1993. Note: Water activity of Spanish intermediate moisture fish products. *Revista Espanola De Ciencia Y Tecnologia De Alimentos* 33:651-656.
- Hand, L. 1994. Controlling water activity and pH in snack sticks. *Meat Marketing and Technology* May:55-56.
- Lee, M.B., and S. Styliadis. 1996. A survey of pH and water activity levels in processed salamis and sausages in Metro Toronto. *Journal of Food Protection* 59:1007-1010.
- Luecke, F.K. 1994. Fermented meat products. *Food Res Intl* 27:299-307.
- Minegishi, Y., Y. Tsukamasa, K. Miake, T. Shimasaki, C. Imai, M. Sugiyama, and H. Shinano. 1995. Water activity and microflora in commercial vacuum-packed smoked salmons. *Journal of the Food Hygienic Society of Japan* 36:442-446.
- Nunez, F., M.C. Diaz, M. Rodriguez, E. Aranda, A. Martin, and M.A. Asensio. 2000. Effects of substrate, water activity, and temperature on growth and verrucosidin production by *Penicillium polonicum* isolated from dry-cured ham. *Journal of Food Protection* 63:231-236.
- Placido, M. and M.P. Aleman. 2002. Rapid hygrometric method for determining water activity. *Ciencia y Tecnologia Alimentaria* 3(4):229-235.

- Rocha-Garza, A.E., and J.F. Zayas. 1996. Quality of broiled beef patties supplemented with wheat germ protein our. *Journal of Food Science* 61:418-421
- Sabadini, E., M.D. Hubinger, P.-J.d.Sobral, and B.C. Carvalho, Jr. 2001. Change of water activity and meat colour in the elaboration process of dehydrated salted meat. *Ciencia e Tecnologia de Alimentos* 21(1):14-19.
- Shimasaki, T., K. Miake, Y. Tsukamasa, M.A. Sugiyama, Y. Minegishi, and H. Shinano. 1994. Effect of water activity and storage temperature on the quality and microora of smoked salmon. *Nippon Suisan Gakkaishi* 60:569-576.
- Untermann, F., and C. Muller. 1992. Influence of aw value and storage temperature on the multiplication and enterotoxin formation of staphylococci in dry-cured raw hams. *International Journal of Food Microbiology* 16:109-115.
- Williams, S.K., G.E. Rodrick, and R.L. West. 1995. Sodium lactate affects shelf life and consumer acceptance of fresh Catfish (*Ictalurus nebulosus*, *marmoratus*) fillets under simulated retail conditions. *Journal of Food Science* 60:636-639.

8.3.2 Mléčné produkty

- Clavero, M.R.S., and L.R. Beuchat. 1996. Survival of *Escherichia coli* O157:H7 in broth and processed salami as influenced by pH, water activity, and temperature and suitability of media for its recovery. *Appl Environ Microbiol* 62:2735-2740.
- Correia, R., M. Magalhaes, M. Pedrini, A. da Cruz, and I. Clementino. 2008. Ice cream made from cow and goat milk: chemical composition and melting point characteristics. *Revista Ciencia Agronomica* 39:251-256.
- Duffy, L.L., P.B. Vanderlinde, and F.H. Grau. 1994. Growth of *Listeria monocytogenes* on vacuum-packed cooked meats: Effects of pH, aw, nitrite and ascorbate. *International Journal of Food Microbiology* 23:377-390.
- Gmez, R., and J. Fernandez-Salguero. 1993. Note: Water activity of Spanish intermediate moisture fish products. *Revista Espanola De Ciencia Y Tecnologia De Alimentos* 33:651-656.
- Hand, L. 1994. Controlling water activity and pH in snack sticks. *Meat Marketing and Technology* May:55-56.
- Hardy, J., J. Scher, and S. Banon. 2002. Water activity and hydration of dairy powders. *Lait* 82:441-442.
- Lee, M.B., and S. Styliadis. 1996. A survey of pH and water activity levels in processed salamis and sausages in Metro Toronto. *Journal of Food Protection* 59:1007-1010.
- Luecke, F.K. 1994. Fermented meat products. *Food Res Intl* 27:299-307.
- Malec, L.S., A.S. Pereyra-Gonzales, G.B. Naranjo, and M.S. Vigo. 2002. Influence of water activity and storage temperature on lysine availability of a milk like system. *Food Res Intl* 35(9):849-853.
- Minegishi, Y., Y. Tsukamasa, K. Miake, T. Shimasaki, C. Imai, M. Sugiyama, and H. Shinano. 1995. Water activity and microflora in commercial vacuum-packed smoked salmons. *Journal of the Food Hygienic Society of Japan* 36:442-446.
- Rocha-Garza, A.E., and J.F. Zayas. 1996. Quality of broiled beef patties supplemented with wheat germ protein our. *Journal of Food Science* 61:418-421.

- Shah, N.P., and R.R. Ravula. 2000. Influence of water activity on fermentation, organic acids production and viability of yoghurt and probiotic bacteria. *Australian Journal of Dairy Technology* 55(3):127-131.
- Shimasaki, T., K. Miake, Y. Tsukamasa, M.A. Sugiyama, Y. Minegishi, and H. Shinano. 1994. Effect of water activity and storage temperature on the quality and microora of smoked salmon. *Nippon Suisan Gakkaishi* 60:569-576.
- Untermann, F., and C. Muller. 1992. Influence of aw value and storage temperature on the multiplication and enterotoxin formation of staphylococci in dry-cured raw hams. *International Journal of Food Microbiology* 16:109-115.
- Williams, S.K., G.E. Rodrick, and R.L. West. 1995. Sodium lactate affects shelf life and consumer acceptance of fresh Catfish (*Ictalurus nebulosus*, *marmoratus*) fillets under simulated retail conditions. *Journal of Food Science* 60:636-639.

8.3.3 Ovoce a zelenina

- Ayub, M., R. Khan, S. Wahab, A. Zeb, and J. Muhammad. 1995. Effect of crystalline sweeteners on the water activity and shelf stability of osmotically dehydrated guava. *Sarhad Journal of Agriculture* 11:755-761.
- Beveridge, T., and S.E. Weintraub. 1995. Effect of blanching pretreatment on color and texture of apple slices at various water activities. *Food Res Intl* 28:83-86.
- Clavero, M.R.S., R.E. Brackett, L.R. Beuchat, and M.P. Doyle. 2000. Influence of water activity and storage conditions on survival and growth of proteolytic *Clostridium botulinum* in peanut spread. *Food Microbiology* 17(1):53-61.
- Fouskaki, M., K. Karametsi, and N.A. Chaniotakis. 2003. Method for the determination of water content in sultana raisins using a water activity probe. *Food Chem* 82:133-1337.
- Gogus, F., C. Cuzdemir, and S. Eren. 2000. Effects of some hydrocolloids and water activity on nonenzymic browning of concentrated orange juice. *Nahrung* 44(6):438-442.
- Hubinger, M., F.C. Menegalli, R.J. Aguerre, and C. Suarez. 1992. Water vapor adsorption isotherms of guava, mango and pineapple. *Journal of Food Science* 57:1405-1407.
- Jimenez, M., M. Manez, and E. Hernandez. 1996. Influence of water activity and temperature on the production of zearalenone in corn by three *Fusarium* species. *International Journal of Food Microbiology* 29:417-421.
- Khalloufi, S., J. Giasson, and C. Ratti. 2000. Water activity of freeze dried mushrooms and berries. *Canadian Agricultural Engineering* 42(1):51-56.
- Kiranoudis, C.T., Z.B. Maroulis, E. Tsami, and D. Marinos-Kouris. 1993. Equilibrium moisture content and heat of desorption of some vegetables. *Journal of Food Engineering* 20:55-74.
- Lopez-Malo, A., and E. Palou. 2000. Modeling the growth/nogrowth interface of *Zygosaccharomyces bailii* in Mango puree. *Journal of Food Science*: 65:516-520.
- Makower, B., and S. Myers. 1943. A new method for the determination of moisture in dehydrated vegetables. *Proceedings of Institute of Food Technologists, 4th Conference* 156.
- Maltini, E., D. Torreggiani, B.R. Brovotto, and G. Bertolo. 1993. Functional properties of reduced moisture fruits as ingredients in food systems. *Food Res Intl* 26:413-419.

- Marin, S., N. Magan, M. Abellana, R. Canela, A.J. Ramos, and V. Sanchis. 2000. Selective Effect of propionates and water activity on maize mycoora and impact on fumonisin B1 accumulation. *Journal of Stored Products Research* 36:203-214.
- Marin, S., V. Sanchis, I. Vinas, R. Canela, and N. Magan. 1995. Effect of water activity and temperature on growth and fumonisin B-1 and B-2 production by *Fusarium proliferatum* and *F. moniliforme* on maize grain. *Lett Appl Microbiol* 21:298-301.
- Monsalve-Gonzalez, A., G.V. Barbosa-Canovas, and R.P. Cavalieri. 1993. Mass transfer and textural changes during processing of apples by combined methods. *Journal of Food Science* 58:1118-1124.
- Pinsirodom, P., and K.L. Parkin. 2000. Selectivity of Celite immobilized patatin (lipid acyl hydrolase) from potato (*Solanum tuberosum* L.) tubers in esterification reactions as influenced by water activity and glycerol analogues as alcohol acceptors. *J. Agric. Food Chem.* 48(2):155-160.
- Tapia de Daza, M.S., C.E. Aguilar, V. Roa, and R.V. Diaz de Tablante. 1995. Combined stress Effects on growth of *Zygosaccharomyces rouxii* from an intermediate moisture papaya product. *Journal of Food Science* 60:356-359.
- Zeb, A., R. Khan, A. Khan, M. Saeed, and S.A. Manan. 1994. Influence of crystalline sucrose and chemical preservatives on the water activity and shelf stability of intermediate banana chips. *Sarhad Journal of Agriculture* 10:721-726.
- Zhang, X.W., X. Liu, D.X. Gu, W. Zhou, R.L. Wang, and P. Liu. 1996. Desorption isotherms of some vegetables. *Journal of the Science of Food and Agriculture* 70:303-306.

8.3.4 Pečivo a cereálie

- Abellana, M., A.J. Ramos, V. Sanchis, and P.V. Nielsen. 2000. Effect of modified atmosphere packaging and water activity on growth of *Eurotium amstelodami*, *E. chevalieri* and *E. herbariorum* on a sponge cake analogue. *Journal of Applied Microbiology* 88:606-616.
- Aramouni, F.M., K.K. Kone, J.A. Craig, and D.Y.C. Fung. 1994. Growth of *Clostridium sporogenes* PA 3679 in home-style canned quick breads. *Journal of Food Protection* 57:882-886.
- Cahagnier, B., L. Lesage, and D. Richard-Molard. 1993. Mould growth and conidiation in cereal grains as affected by water activity and temperature. *Lett Appl Microbiol* 17:7-13.
- Clawson, A.R., and A.J. Taylor. 1993. Chemical changes during cooking of wheat. *Food Chem* 47:337-343.
- Fleurat-Lessard, F. 2002. Qualitative reasoning and integrated management of the quality of stored grain: a promising new approach. *Journal of Stored Products Research* 38:191-218.
- Gmez, R., J. Fernandez-Salguero, M.A. Carmona, and D. Sanchez. 1993. Water activity in foods with intermediate moisture levels: Bakery and confectionery products: Miscellany. *Alimentaria* 30:55-57.
- Guynot, M.E., A.J. Ramos, L. Seto, P. Purroy, V. Sanchis, and S. Marin. 2003. Antifungal activity of volatile compounds generated by essential oils against fungi commonly causing deterioration of bakery products.
- Harris, M., and M. Peleg. 1996. Patterns of textural changes in brittle cellular cereal foods caused by moisture sorption. *Cereal Chem* 73:225-231.

- Hope, R., and N. Magan. 2003. Two-dimensional environmental profiles of growth, deoxynivalenol and nivalenol production by *Fusarium culmorum* on wheat-based substrate. *Lett Appl Microbiol* 37:70-74.
- Michniewicz, J., C.G. Biliaderis, and W. Bushuk. 1992. Effect of added pentosans on some properties of wheat bread. *Food Chem* 43:251-257.
- Moreno-Contreras, M.D., A.J. Martinez-Yepep, and R.R. Martinez. 2000. Determination of deoxynivalenol (DON) in wheat, barley and corn and its relationship with the levels of total molds, *Fusarium* spp., infestation percentage, and water activity. *Archivos Latinoamericanos de Nutricion*. 50(2):183-186.
- Phoungchandang, S., and J.L. Woods. 2000. Moisture diffusion and desorption isotherms for banana. *Journal of Food Science* 65:651-657.
- Ramanathan, S., and S. Cenkowski. 1995. Sorption isotherms of flour and flow behaviour of dough as influenced by our compaction. *Canadian Agricultural Engineering* 37:119-124.
- Roessler, P.F., and M.C. Ballenger. 1996. Contamination of an unpreserved semisoft baked cookie with a Xerophilic *Aspergillus* species. *Journal of Food Protection* 59:1055-1060.
- Schebor, C., and J. Chirife. 2000. A survey of water activity and pH values in fresh pasta packed under modified atmosphere manufactured in Argentina and Uruguay. *Journal of Food Protection* 63:965-969.
- Seiler, D.A.L. 1979. The mould-free shelf life of bakery products. *FMBRA Bulletin* April:71-74.
- Sumner, S.S., J.A. Albrecht, and D.L. Peters. 1993. Occurrence of enterotoxigenic strains of *Staphylococcus aureus* and enterotoxin production in bakery products. *Journal of Food Protection* 56:722-724.
- Tesch, R., M.D. Normand, and M. Peleg. 1996. Comparison of the acoustic and mechanical signatures of two cellular crunchy cereal foods at various water activity levels. *Journal of the Science of Food and Agriculture* 70:347-354.
- Weegels, P.L., J.A. Verhoek, A.M.G. de Groot, and R.J. Hamer. 1994. Effects of gluten of heating at different moisture contents: I. Changes in functional properties. *Journal of Cereal Science* 19:31-38.

8.3.5 Nápoje, polévky, omáčky a konzervy

- Cardelli, C., and T.P. Labuza. 2001. Application of Weibull Hazard Analysis to the determination of shelf life of roasted and ground coffee. *Lebensm Wiss Technol* 34:273-278.
- Carson, K.J., J.L. Collins, and M.P. Penfield. 1994. Unrefined, dried apple pomace as a potential food ingredient. *Journal of Food Science* 59:1213-1215.
- Cavia, M.M., M.A. Fernandez-Muio, J.F. Huidobro, and M.T. Sancho. 2004. Correlation between Moisture and Water Activity of Honeys Harvested in Different Years. *Journal of Food Science* 69:C-368-370.
- Durrani, M.J., R. Khan, M. Saeed, and A. Khan. 1992. Development of concentrated beverages from Anna apples with or without added preservatives by controlling activity of water for shelf stability. *Sarhad Journal of Agriculture* 8:23-28.
- Ferragut, V., J.A. Salazar, and A. Chiralt. 1993. Stability in the conservation of emulsified sauces low in oil content. *Alimentaria* 30:67-69.
- Gleiter, R.A., H. Horn, and H.-D. Isengard. 2006. Influence of type and state of crystallization on the water activity of honey. *Food Chem* 96:441-445.

- Hajmeer, M.N., F.M. Aramouni, and E.A.E.Boyle. 2000. Shelflife of lite syrup after opening and storage at room or refrigerated temperature. *Journal of Food Quality* 23:529-540.
- Ibarz, A., J. Pagan, and R. Miguelsanz. 1992. Rheology of clarified fruit juices: II. Blackcurrant juices. *Journal of Food Engineering* 15:63-74.
- Khalloufi, S., Y. El-Maslouhi, and C. Ratti. 2000. Mathematical model for prediction of glass transition temperature of fruit powders. *Journal of Food Science* 65:842-848.
- Kusumegi,K., T.Takahashi, and M.Miyagi. 1996. Effects of addition of sodium citrate on the pasteurizing conditions in Tuyu, Japanese noodle soup. *Journal of the Japanese Society for Food Science and Technology* 43:740-747.
- Perera, C.O. 2005. Selected quality attributes of dried foods. *Drying Technology* 23:717-730.
- Sa, M.M., and A.M. Sereno. 1993. Effect of temperature on sorption isotherms and heats of sorption of quince jam. *International Journal of Food Science & Technology* 28:241-248.
- Shafi ur-Rahman, M. 2005. Dried food properties: challenges ahead. *Drying Technology* 23:695-715.

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- Ahlneck, C., and G. Zografi. 1990. The Molecular basis of moisture effects on the physical and chemical stability of drugs in the solid state. *International Journal of Pharmaceutics* 62:87-95.
- Bell, L.N., and K.L. White. 2000. Thiamin Stability in Solids as Affected by the Glass Transition. *Journal of Food Science* 65:498-501.
- Cochet, N., and A.L. Demain. 1996. Effect of water activity on production of beta-lactam antibiotics by *Streptomyces clavuligerus* in submerged culture. *Journal of Applied Bacteriology* 80:333-337.
- Constantino, H.R., R. Langer, and A.M. Klibanov. 1994. Solid-Phase Aggregation of Proteins under Pharmaceutically Relevant Conditions. *Journal of Pharmaceutical Science* 83:1662-1669.
- Enigl, D.C. 2001. Pharmaceutical stability testing using water activity. *European Pharmaceutical Review* 6:46-49.
- Enigl, D.C., and K.M.Sorrel. 1997. Water Activity and Self-Preserving Formulas. p. 45-73. In J.J. Kabara, and D.S. Orth (ed.) *Preservative-Free and Self-Preserving Cosmetics and Drugs: Principles and Practice*. Marcel Dekker.
- Hageman, M.J. 1988. The Role of Moisture in Protein Stability. *Drug Dev. Ind. Pharm.* 14:2047-2070.
- Heidemann, D.R., and P.J. Jarosz. 1991. Preformulation Studies Involving Moisture Uptake in Solid Dosage Forms. *Pharmaceutical Research* 8:292-297.
- Kontny, M.J. 1988. Distribution of Water in Solid Pharmaceutical Systems. *Drug Dev. Ind. Pharm.* 14:1991-2027.
- Sablani, S.S., K. Al-Belushi, I. Al-Marhubi, and R. Al-Belushi. 2007. Evaluating Stability of Vitamin C in Fortified Formula Using Water Activity and Glass Transition. *International Journal of Food Properties* 10:61-71.
- Zografi, G. 1988. States of Water Associated with Solids. *Drug Dev. Ind. Pharm.* 14:1905-1926.

Zografi, G., and M.J. Kontny. 1986. The interactions of water with cellulose- and starch-derived pharmaceutical excipients. *Pharmaceutical Research* 3:187-193.

8.5 Různé

- Bell, L.N. 1995. Kinetics of non-enzymatic browning in amorphous solid systems: Distinguishing the Effects of water activity and the glass transition. *Food Res Intl* 28:591-597.
- Bell, L.N., and T.P. Labuza. 1992. Compositional influence on the pH of reduced-moisture solutions. *Journal of Food Science* 57:732-734.
- Bell, L.N., and T.P. Labuza. 1994. Influence of the low-moisture state on pH and its implication for reaction kinetics. *Journal of Food Engineering* 22:291-312.
- Bhandari, B., and I. Bareyre, 2003. Estimation of crystalline phase present in glucose crystal-solution mixture by water activity measurement. *Lebensm Wiss Technol* 36:729-733(5).
- Brake, N.C., and O.R. Fennema. 1993. Edible coatings to inhibit lipid migration in a confectionery product. *Journal of Food Science* 58:1422-1425.
- Dole, M., and L. Faller. 1950. Water sorption by synthetic high polymers. *Journal of the American Chemical Society* 72:414-419.
- Fernandez-Salguero, J., R. Gmez, and M.A. Carmona. 1993. Water activity in selected high-moisture foods. *Journal of Food Composition and Analysis* 6:364-369.
- Juhan, K., and G.K. Byung. 2000. Lipase-catalyzed synthesis of lysophosphatidylcholine using organic cosolvent for in situ water activity control. *Journal of American Oil Chemists' Society* 77(7):701-797.
- Lima, J.R., S.D.S. Campos, and L.-A.G. Goncalves. 2000. Relationship between water activity and texture of roasted and salted cashew kernel. *Journal of Food Science and Technology* 37(5):512-513.
- Lomauro, C.J., A.S. Bakshi, and T.P.Labuza. 1985a. Evaluation of food moisture sorption isotherm equations. Part II: Milk, coffee, tea, nuts, oilseeds, spices and starchy foods. *Lebensm Wiss Technol* 18:118-124.
- Lomauro, C.J., A.S. Bakshi, and T.P. Labuza. 1985b. Evaluation of food moisture sorption isotherm equations. Part I: Fruit, vegetable and meat products. *Lebensm Wiss Technol* 18:111-117.

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Copyright ©2007-2016
All rights reserved